MAS® DOA-XSE®

LIQUID ASSAYED DRUGS OF ABUSE CONTROL

For In Vitro Diagnostic Use

INTENDED USE

Multi-analyte **MAS® DOA-XSE®** is intended for use as a consistent test sample for monitoring assay conditions in semi-quantitative and qualitative analyses of patient urine specimens for drugs and drug metabolites. Include the controls with patient urine specimens when assaying for any of the listed constituents.

PRODUCT DESCRIPTION

DOA-XSE® is prepared from certified drug free human urine pools. Analyte levels are adjusted with purified drugs or drug metabolites. Preservatives and stabilizers are added to maintain product integrity.

Level 1 is formulated to assure a negative analytical response for all assay methods (i.e. below the level of quantitation). **Level 2** and **Level 3** bracket immunoassay procedure cutoffs. The cutoff levels for Amphetamines, PCP (Phencyclidine), Opiates, Cocaine and Marijuana (Cannabinoid) are applied in accordance with the Mandatory Guidelines for Federal Workplace Drug Testing programs /1/ as established by the Substance Abuse and Mental Health Services Administration (SAMHSA). **Level 4** is an elevated positive control. Although the constituent levels in **Level 2** and **Level 3** are targeted for specific SAMHSA cutoffs as appropriate, the controls may be used with all semi-quantitative or qualitative procedures where the levels are appropriate for the user's needs.

Good laboratory practice dictates that each laboratory establish its own acceptable values for these controls. Recovered values are dependent upon the method used, standardization, method interferences and/or cross reactivities and antibody specificity. Refer to reagent kit manufacturer's recovery claims as an aid in establishing values for specific methods used.

CAUTION: DOA-XSE® is prepared from human urine pools. Each pool has been tested using FDA accepted methods and found to be non-reactive for Hepatitis B surface antigen (HBsAg), Hepatitis C Antibody (HCV), and HIV-1 and HIV-2 antibody. However, no test method can offer complete assurance that products derived from human source material are free of infectious agents. This control must be handled in accordance with recommendations from Centers for Disease Control/National Institutes of Health manual, "Biosafety in Microbiological and Biomedical Laboratories," 1999.

DOA-XSE® contains small quantities of sodium azide (<0.1%) which may react with lead or copper plumbing to form metal azides. Continuous buildup of metal azides can present an explosion hazard. Upon disposal, flush drains with copious amounts of water to prevent azide buildup.

CAUTION: The packaging of this product contains dry natural rubber.

STORAGE AND STABILITY

Store **DOA-XSE®** at 2-8°C. Unopened vials are stable until the expiration date on the label. Once opened, vials of control are stable for 30 days when stored tightly capped at 2-8°C. **Do not freeze.**

Bacterial contamination produces an increase in turbidity and/or a characteristic odor. Discard vial if evidence of microbial contamination is observed.

INSTRUCTIONS FOR USE

Thoroughly mix the contents of the vial before each use by gently inverting for several minutes. Open the vial and transfer the required quantity of control into a clean sample cup. Replace the cap immediately and store the opened vial at 2-8°C. Assay controls in accordance with the reagent and/or instrument manufacturer's instructions for unknown specimens.

Include **DOA-XSE®** in the test run following the same procedures used for patient samples.

LIMITATIONS OF PROCEDURE

It has been reported that certain types of plastic and glass laboratory ware will bind the cannabinoid metabolite present in these controls. Use new pipettes and sample cups for each patient sample, calibrator and control. Assay samples without delay.

The concentration of certain constituents in **Level 4** may exceed the assay range of some analytical methods. Dilute **Level 4**, as necessary, with **Level 1** (Negative) **DOA-XSE®**. Results obtained must be corrected for dilution.

Variability in antibody specificity with different reagent lots may affect values obtained with immunological procedures.

Accurate and reproducible results are dependent upon properly functioning instruments, reagents, and good laboratory technique. This product is intended for use as an assayed control for quantitative assays of listed constituents in human serum. This product is not intended for use as a calibrator.

TECHNICAL ASSISTANCE

In the USA, for technical assistance, insert updates and information, call 800-232-3342 or 510-979-5417. Outside the USA, please contact your local authorized distributor.

BIBLIOGRAPHY

/1/ Substance Abuse and Mental Health Services Administration, Mandatory Guidelines for Federal Workplace Drug Testing Programs, Federal Register, June 9, 1994, 59 FR, 29916.

Cat. No.	Description	Size
SDOA-101	Level 1, Negative	6 x 18 mL
SDOA-202	Level 2, Low challenge to immunoassay cutoff	6 x 18 mL
SDOA-303	Level 3, High challenge to immunoassay cutoff	6 x 18 mL
SDOA-404	Level 4, Elevated immunoassay	6 x 18 mL

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(a) Weighed-in amounts are not relevant due to inherent process loss.

- (a) Las cantidades pesadas no son pertinentes debido a la pérdida inherente al proceso.
- (a) Eingewogene Mengen wegen natürlichen Verfahrensverlustes irrelevant.
- (a) Les quantités pesées ne sont pas significatives du fait de la perte inhérente au procédé.
- (a) I valori ponderati non sono rilevanti a causa di perdite intrinseche del processo.

(b) No GC/MS method available for Ethanol quantitation. Confirmation values for Ethanol obtained using Hitachi.

- (b) No existe un método por GC/MS para cuantificar el etanol. Los valores de confirmación del etanol se han obtenido usando Hitachi.
- (b) Keine GC/MS-Methode zur Ethanolmengenbestimmung erhältlich. Bestätigungswerte für Ethanol mit Hitachi festgelegt.
- (b) Pas de méthode GC/MS disponible pour la quantification de l'éthanol. Les valeurs de confirmation concernant l'éthanol ont été obtenues par méthode Hitachi.
- (b) Non è disponibile alcun metodo di gas cromatografia/spettrometria di massa per la quantificazione dell'etanolo. I valori di conferma per l'etanolo sono stati ottenuti con Hitachi.

**** Theoretical target for Hitachi®**


- ** Objetivo teórico para Hitachi®
- ** Theoretisches Ziel für Hitachi®
- ** Cible théorique pour Hitachi®
- ** Obiettivo teorico per Hitachi®

Hitachi® Reg TM: Roche Diagnostics, Indianapolis, IN

A comprehensive range of Quality Controls for your Laboratory/ Un rango comprensivo de control de calidad para su laboratorio /Eine umfassende Palette von Qualitätskontrollprodukten für Ihr Labor/ Un menu complet de contrôles de qualité pour votre laboratoire/ Range globale degli analiti per Controllo di Qualità del vostro Laboratorio

Serum Chemistry/ La Química del Suero/ Serumchemie/ Chimie de Sérum/ Chimica di Serum	
MAS[®] chemTRAK[®] · H LIQUID HUMAN BASED CHEMISTRY CONTROL/ CONTROL LÍQUIDO PARA BIOQUÍMICA/ CHEMIE-KONTROLLFLÜSSIGKEIT/ SOLUTION DE CONTROLE DE QUALITE LIQUIDE/ CONTROLLO LIQUIDO PER LA CHIMICA	MAS[®] Bilirubin LIQUID ASSAYED BILIRUBIN CONTROL/ CONTROL LÍQUIDO VALORADO PARA BILIRUBINA/ ANALYSIERTE BILIRUBIN-KONTROLLFLÜSSIGKEIT/ CONTROLE DE QUALITE LIQUIDE POUR LE DOSAGE DE BILIRUBINE/ CONTROLLO LIQUIDO TESTATO PER BILIRUBINA
MAS[®] chemTRAK[®] · Platinum LIQUID CHEMISTRY CONTROL/ CONTROL LÍQUIDO PARA BIOQUÍMICA/ CHEMIE-KONTROLLFLÜSSIGKEIT/ SOLUTION DE CONTROLE DE QUALITE LIQUIDE/ CONTROLLO LIQUIDO PER LA CHIMICA	MAS[®] Alcohol/Ammonia LIQUID ASSAYED ALCOHOL/AMMONIA CONTROL/ CONTROL LÍQUIDO VALORADO PARA ALCOHOL/AMONÍACO/ ANALYSIERTE ALCOHOL/AMMONIA-KONTROLLFLÜSSIGKEIT/ CONTROLE DE QUALITE LIQUIDE POUR LE DOSAGE D'ALCOOL ET D'AMMONIAQUE/ CONTROLLO LIQUIDO TESTATO PER ALCOOL/AMMONIACA
MAS[®] Lipid LIQUID ASSAYED LIPID CONTROL/ CONTROL LÍQUIDO ANALIZADO DE LÍPIDOS/ FLUSSIGE ANALYSIERTE LIPIDKONTROLLE/ CONTROLE LIQUIDE POUR DOSAGE DES LIPIDES/ CONTROLLO CAMPIONATO LIQUIDO LIPIDI	MAS[®] CSF LIQUID ASSAYED CEREBROSPINAL FLUID CONTROL/ CONTROL LÍQUIDO VALORADO PARA LÍQUIDO CEFALORRAQUÍDEO/ ANALYSIERTE LIQUOR CEREBROSPINALIS-KONTROLLFLÜSSIGKEIT/ CONTROLE DE QUALITE LIQUIDE POUR LES DOSAGES DANS LE L.C.R./ CONTROLLO LIQUIDO TESTATO PER LIQUIDO CEREBROSPINALE
Urine Chemistry/ Bioquímica de la Orina/ Urinchemie/ Dosage des Urines/ Chimica delle Urine	
MAS[®] UrichemTRAK LIQUID ASSAYED URINE CONTROL/ CONTROL LÍQUIDO VALORADO PARA ORINA/ ANALYSIERTE URIN-KONTROLLFLÜSSIGKEIT/ CONTROLE DE QUALITE LIQUIDE POUR LE DOSAGE DES URINES/ CONTROLLO LIQUIDO TESTATO PER LE URINE	MAS[®] UA Control LIQUID ASSAYED URINALYSIS CONTROL/ CONTROL LÍQUIDO VALORADO PARA ORINA/ ANALYSIERTE URIN-KONTROLLFLÜSSIGKEIT/ CONTROLE DE QUALITE LIQUIDE POUR LE DOSAGE DES URINALYSIS/ CONTROLLO LIQUIDO TESTATO PER LE URINALYSIS
Immunoassay/ Inmunoanálisis/ Immunoassay/ Immunologique/ Immunologici	
MAS[®] Liquimmune[®] LIQUID ASSAYED IMMUNOASSAY CONTROL/ CONTROL LÍQUIDO VALORADO PARA INMUNOANÁLISIS/ ANALYSIERTE IMMUNOASSAY-KONTROLLFLÜSSIGKEIT/ CONTROLE DE QUALITE LIQUIDE POUR LE DOSAGE IMMUNOLOGIQUE/ CONTROLLO LIQUIDO TESTATO PER SAGGI IMMUNOLOGICI	MAS[®] T-Marker LIQUID ASSAYED TUMOR MARKER CONTROL/ CONTROL LÍQUIDO VALORADO PARA TUMOR MARKER/ ANALYSIERTE TUMOR MARKER-KONTROLLFLÜSSIGKEIT/ CONTROLE DE QUALITE LIQUIDE POUR LE MARQUEUR TUMORAUX/ CONTROLLO LIQUIDO TESTATO PER MARCATORI TUMORALI
MAS[®] PAR TDM LIQUID ASSAYED TDM CONTROL/ CONTROL LÍQUIDO VALORADO PARA FÁRMACOS/ ANALYSIERTE THERAPIEMEDIKATION-KONTROLLFLÜSSIGKEIT/ CONTROLE DE QUALITE LIQUIDE POUR LE DOSAGE DES MEDICAMENTS THERAPEUTIQUES/ CONTROLLO LIQUIDO TESTATO PER FARMACI TERAPEUTICI	
Cardiac/ Cardíacos/ Herz/ Cardiaques/ Cardiaci	
MAS[®] Cardioimmune[®] - XL LIQUID ASSAYED CARDIAC MARKER CONTROL/ CONTROL LÍQUIDO VALORADO PARA MARCADORES CARDÍACOS/ ANALYSIERTE HERZMARKER-KONTROLLFLÜSSIGKEIT/ CONTROLE DE QUALITE LIQUIDE POUR LE DOSAGE DE MARQUEURS CARDIAQUES/ CONTROLLO LIQUIDO TESTATO PER MARCATORI CARDIACI	
Protein/ Proteína/ Eiweiß/ Protéine/ Proteina	
MAS[®] Immunology LIQUID ASSAYED IMMUNOLOGY CONTROL/ CONTROL LÍQUIDO VALORADO PARA IMMUNOLOGÍA/ ANALYSIERTE IMMUNOLOGY-KONTROLLFLÜSSIGKEIT/ CONTROLE DE QUALITE LIQUIDE POUR LE DOSAGE D'IMMUNOLOGIE/ CONTROLLO LIQUIDO TESTATO PER IMMUNOLOGIA	
Specialty/ Especialidad/ Spezialität/ Spécialité/ Specialità	
MAS[®] Diabetes LIQUID ASSAYED DIABETES CONTROL/ CONTROL LÍQUIDO VALORADO PARA DIABETES / ANALYSIERTE DIABETES-KONTROLLFLÜSSIGKEIT / CONTROLE DE QUALITE LIQUIDE POUR LE DIABETE /CONTROLLO LIQUIDO TESTATO PER IL DIABETE	
Toxicology/ Toxicología/ Toxikologie/ Toxicologie/ Toxicology	
MAS[®] DOA-XSE[®] • MAS[®] DOA-XAB[®] • MAS[®] DOA GC/MS • MAS[®] Opiates 2000 • MAS[®] Opiates 2000 GC/MS LIQUID ASSAYED DRUGS OF ABUSE CONTROLS/ CONTROL LÍQUIDO VALORADO PARA DROGAS DE ABUSO/ ANALYSIERTE DROGENMISSBRAUCH-KONTROLLFLÜSSIGKEIT/ CONTROLE DE QUALITE LIQUIDE POUR LE DOSAGE DROGUES D'ABUS/ CONTROLLO LIQUIDO TESTATO PER DROGHE D'ABUSO	

Rev. 2007-1

LOT	Lot Number / Número de Lote / Chargennummer / Numéro de lot / Numero di lotto
	"Use By" date / Fecha de caducidad / Haltbarkeitsdatum / Date limite d'utilisation optimale / "Utilizzare entro"
	Refrigerated temperature limitation / Limitación de la temperatura de refrigeración / Kühltemperaturgrenze / Limites de température pour un colis réfrigéré / Limite delle temperature di raffreddamento
REF	Catalog Number / Número De Catálogo / Katalognummer / Catalogue N° / Numero di catalogo
	Manufacturer / Fabricante / Hersteller / Fabricant / Produttore
	Caution / Precaución / Vorsicht / Attention / Attenzione
	Consult instructions for use / Consultar las instrucciones de uso / Gebrauchsanweisung beachten / Se référer au mode d'emploi / Consultare le istruzioni per l'uso
EC REP	Authorized Representative / Representante autorizado / Vertragshändler / Mandataire / Rappresentante autorizzato
IVD	For In Vitro Diagnostic Use / Para uso en diagnóstico in vitro / In-vitro-Diagnostikum / Diagnostics in vitro / Per uso diagnostico in vitro

EC REP

Microgenics GmbH
Spitalhofstrasse 94
D-94032 Passau
Germany
Tel: +49 (0) 851-88 68 90
Fax: +49 (0) 851-88 68 910

Manufacturer / Fabricante / Hersteller / Fabricant / Produttore:
MICROGENICS CORPORATION
46360 Fremont Blvd.
Fremont, CA 94538 USA
www.microgenics.com
ISO 13485 Certified Company

microgenics **MAS**
CONTROLS

MAS[®] DOA-XSE[®]
Liquid Assayed Drugs of abuse Control

6 x 18mL

IVD

REF SDOA-101

LOT XSE09041

2009-04-30

Microgenics Corporation
46360 Fremont Blvd.
Fremont, CA 94538 USA

CE

EC REP

Microgenics GmbH
Spitalhofstrasse 94
D-94032 Passau
Germany

XSE-INS-1
Rev. 6

8
XSE-INS
Rev. 6

microgenics **MAS**
CONTROLS

MAS[®] DOA-XSE[®]
Liquid Assayed Drugs of abuse Control

4 x 18mL

TK 59-27

IVD

REF SDOA-S1

LOT XSE0904S

2009-04-30

Microgenics Corporation
46360 Fremont Blvd.
Fremont, CA 94538 USA

CE

EC REP

Microgenics GmbH
Spitalhofstrasse 94
D-94032 Passau
Germany

XSE-INS-S
Rev. 6

microgenics **MAS**
CONTROLS

MAS[®] DOA-XSE[®]
Liquid Assayed Drugs of abuse Control

6 x 18mL

TK 59-24

IVD

REF SDOA-202

LOT XSE09042

2009-04-30

Microgenics Corporation
46360 Fremont Blvd.
Fremont, CA 94538 USA

CE

EC REP

Microgenics GmbH
Spitalhofstrasse 94
D-94032 Passau
Germany

XSE-INS-2
Rev. 6

microgenics **MAS**
CONTROLS

MAS[®] DOA-XSE[®]
Liquid Assayed Drugs of abuse Control

6 x 18mL

TK 59-25

IVD

REF SDOA-303

LOT XSE09043

2009-04-30

Microgenics Corporation
46360 Fremont Blvd.
Fremont, CA 94538 USA

CE

EC REP

Microgenics GmbH
Spitalhofstrasse 94
D-94032 Passau
Germany

XSE-INS-3
Rev. 6

microgenics **MAS**
CONTROLS

MAS[®] DOA-XSE[®]
Liquid Assayed Drugs of abuse Control

6 x 18mL

TK 59-26

IVD

REF SDOA-404

LOT XSE09044

2009-04-30

Microgenics Corporation
46360 Fremont Blvd.
Fremont, CA 94538 USA

CE

EC REP

Microgenics GmbH
Spitalhofstrasse 94
D-94032 Passau
Germany

XSE-INS-4
Rev. 6